

E. Lopez St., Jaro, Iloilo City "PhilHealth Accredited Health Care Provider"

Tel No.: (033) 320 2431 | Fax No.: (033) 3202623 | Email Address: medcenter@wvsu.edu.ph

SUPPLEMENTAL/BID BULLETIN

Project Reference No. _______ IB No. 24-04 Name of the Project

Procurement of Various Pulmonary Supplies

Location of the Project WVSU Medical Center

ADDENDUM NO. <u>01 - 2023</u>

November 24, 2023

This Addendum No. 01 is issued to modify or amend items in the Bidding Documents. Attached are the Revised Section VII. Technical Specifications (pages 18-21). These shall form an integral part of the Bidding Documents.

A. Section VII. Technical Specifications

ORIGINAL AS INDICATED IN THE BIDDING DOC			ITS	Corrections/Modifications/Amendments (in bold & italicized words)
tem No.	Item Description	Qty	Unit	Item Description
1	BIPAP FULL FACE MASK (NON- VENTED)- MEDIUM	2	PCS	BIPAP FULL FACE MASK (NON- VENTED)- MEDIUM
	• full face mask, vented with anti asphyxiation valve			
	has a tube fit design that maintan an effective seal			
	around nasogastric or orogastric tubes			
	single patient use up to 14 days of use			
	 material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam 			
	 has a simple color coded design to help choose the right mask size 			
	has a dual swivel elbow			
	 can withstand pressures up to 25 cmH20 			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
2	BIPAP FULL FACE MASK (NON- VENTED)- SMALL	3	PCS	BIPAP FULL FACE MASK (NON- VENTED)- SMALL
	• full face mask, vented with anti asphyxiation valve			
	 has a tube fit design that maintan an effective seal Around nasogastric or orogastric tubes 			
	• single patient use up to 14 days of use			
	 material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam 			
	 has a simple color coded design to help choose the right mask size 			
	has a dual swivel elbow			
	 can withstand pressures up to 25 cmH20 			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
3	BIPAP FULL FACE MASK (VENTED)- LARGE	1	PCS	BIPAP FULL FACE MASK (VENTED)- LARGE
	• full face mask, vented with anti asphyxiation valve			
	 has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes 			
	• single patient use up to 14 days of use			2-
	 material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam 			
	 has a simple color coded design to help choose the right mask size 			
	has a dual swivel elbow			
	• can withstand pressures up to 25 cmH20			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation



E. Lopez St., Jaro, Iloilo City "PhilHealth Accredited Health Care Provider"



Tel No.: (033) 320 2431 | Fax No.: (033) 3202623 | Email Address: medcenter@wvsu.edu.ph

SUPPLEMENTAL/BID BULLETIN

Project Reference No. IB No. 24-04 Name of the Project

Procurement of Various Pulmonary Supplies

Location of the Project WVSU Medical Center

ADDENDUM NO. 01 - 2023

November 24, 2023

OR	IGINAL AS INDICATED IN THE BIDDING DO	CUMEN	NTS	Corrections/Modifications/Amendment (in bold & italicized words)		
tem No.	Item Description	Qty	Unit	Item Description		
4	BIPAP FULL FACE MASK (VENTED)- MEDIUM	38	PCS	BIPAP FULL FACE MASK (VENTED)- MEDIUM		
	• full face mask, vented with anti asphyxiation valve					
	 has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes 		-			
	• single patient use up to 14 days of use					
	 material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam 					
	 has a simple color coded design to help choose the right mask size 					
	has a dual swivel elbow					
	• can withstand pressures up to 25 cmH20					
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation		
5	BIPAP FULL FACE MASK (VENTED)- SMALL	38	PCS	BIPAP FULL FACE MASK (VENTED)- SMALL		
	• full face mask, vented with anti asphyxiation valve					
	 has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes 					
	• single patient use up to 14 days of use					
	material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam					
	 has a simple color coded design to help choose the right mask size 					
	has a dual swivel elbow					
	can withstand pressures up to 25 cmH20		-			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation		
6	HEAT AND MOISTURE EXCHANGER (ADULT)	611	PCS	HEAT AND MOISTURE EXCHANGER (ADULT)		
	with mechanical Hydrophobic mebrane			with mechanical Hydrophobic membrane		
	with hygroscopic membrane					
	Individually packaged					
	designed for tidal volume 150-1500 ml			designed for tidal volume 150-1500 ml		
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation		
7	HEAT AND MOISTURE EXCHANGER (PEDIA)	300	PCS	HEAT AND MOISTURE EXCHANGER (PEDIA)		
	with mechanical Hydrophobic mebrane			with mechanical Hydrophobic membrane		
	with hygroscopic membrane					
	Individually packaged					
	designed for tidal volume 75-500ml					
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation		



E. Lopez St., Jaro, Iloilo City "PhilHealth Accredited Health Care Provider"

Tel No.: (033) 320 2431 | Fax No.: (033) 3202623 | Email Address: medcenter@wvsu.edu.ph



SUPPLEMENTAL/BID BULLETIN

Name of the Project

Procurement of Various Pulmonary Supplies

Location of the Project WVSU Medical Center

ADDENDUM NO. 01 - 2023

November 24, 2023

ORIGINAL AS INDICATED IN THE BIDDING DOC			ITS	Corrections/Modifications/Amendments (in bold & italicized words)
tem No.	Item Description	Qty	Unit	Item Description
8	METERED DOSE INHALER PORT	127	PCS	METERED DOSE INHALER PORT
	 diameter- 22 mm (ID), 22 mm (OD) 			
	 with arrow that clearly marks the direction of aiflow 			
	disposable			
	use in the delivery of aerosolized medication			
	individually packed			
	clear and transparent plastic			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
9	PULSE OXIMETER SENSOR, FINGER PROBE (ADULT)	80	PCS	PULSE OXIMETER SENSOR, FINGER PROBE (ADULT)
	disposable			(
	• for single use			110
	 Intended for non-invasive, continuous SpO2 monitoring 			
	individually packed			
	 Ideal for application on the index finger, thumb, or 			
	great toe			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
10	PULSE OXIMETER SENSOR, FINGER PROBE (NEONATE)	80	PCS	PULSE OXIMETER SENSOR, FINGER PROBE (NEONATE)
	disposable		l.	
	for single use			
	individually packed			
	 Intended for non-invasive, continuous SpO2 monitoring 			
	 Ideal for application on the index finger, thumb, or great toe 			
9	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
11	PULSE OXIMETER SENSOR, FINGER PROBE (PEDIA)	80	PCS	PULSE OXIMETER SENSOR, FINGER PROBE (PEDIA)
	disposable	-		
	• for single use	1		
	individually packed	-		
	Intended for non-invasive, continuous SpO2 monitoring			
	 Ideal for application on the index finger, thumb, or great toe 			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation
12	VIRAL FILTER	750	PCS	VIRAL FILTER
	Individually packaged			
	designed for tidal volume 150-1200ml			
	• viral efficiency 99.99+%			
	• 22 mm ID x 15 mm ID/22 mm OD			
	Bacterial efficiency 99.999+%			
	Note: Please provide sample for evaluation			Note: Please provide 1 sample for evaluation

E. Lopez St., Jaro, Iloilo City "PhilHealth Accredited Health Care Provider"



Tel No.: (033) 320 2431 | Fax No.: (033) 3202623 | Email Address: medcenter@wvsu.edu.ph

SUPPLEMENTAL/BID BULLETIN

Project Reference No. _______ IB No. 24-04

Name of the Project

Procurement of Various Pulmonary Supplies

Location of the Project WVSU Medical Center

ADDENDUM NO. 01 - 2023

November 24, 2023

ORIGINAL AS INDICATED IN THE BIDDING DOCUMENTS			Corrections/Modifications/Amendments (in bold & italicized words)		
Item No.	Item Description	Qty	Unit	Item Description	
13	VOLUME RESPIRATOR TUBINGS (ADULT) 500	PCS	VOLUME RESPIRATOR TUBINGS (ADULT)		
	disposable				
	for single use				
	with nebulization kit				
	with water trap				
	dual limb				
	Note: Please provide sample for evaluation	ation		Note: Please provide 1 sample for evaluation	
14	VOLUME RESPIRATOR TUBINGS (PED	DIA) 250	PCS	VOLUME RESPIRATOR TUBINGS (PEDIA)	
	disposable				
	for single use				
	with nebulization kit				
	with water trap				
	dual limb				
	Note: Please provide sample for evaluation	ation		Note: Please provide 1 sample for evaluation	
15	VOLUME RESPIRATOR TUBINGS (NEC	ONATE) 300	PCS	VOLUME RESPIRATOR TUBINGS (NEONATE)	
	disposable				
	for single use				
	 with nebulization kit 				
	with water trap				
	dual limb				
	Note: Please provide sample for evalu	ation		Note: Please provide 1 sample for evaluation	

For guidance and information of all concerned.

Bids and Awards Committee
For the Chairperson:
PATRICIA P. TUPAS, DDM Vice Chairperson Bids and Awards Committee
Received by the Bidder:
Printed Name & Signature
Name of Bidder/Company
Date:

CESAR Q. TRONCO, R.N.

Chairperson

Flocurement of various Pulmonary Supplies Div Unit Brand	Item No.	Specification	Statement of Compliance [Bidders must state here either "Comply" or "Not Comply" agains each of the individual parameters of each Specification stating the corresponding performance parameter of the equipment offered Statements of "Comply" or "Not Comply" must be supported by evidence in a Bidders Bid and cross-referenced to that evidence Evidence shall be in the form of manufacturer's un-amended sales literature, unconditional statements of specification and compliance issued by the manufacturer, samples, independent test data etc., as appropriate. A statement that is not supported by evidence or is subsequently found to be contradicted by the evidence presented will render the Bid under evaluation liable for rejection. A statement either in the Bidder's statement of compliance or the supporting evidence that is found to be false either during Bid evaluation, post- qualification or the execution of the Contract may be regarded as fraudulent and render the Bidder or supplier liable for prosecution subject to the applicable laws and issuances.]				
• full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON-VENTED)- SMALL • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE 1 PCS • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a stube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		Procurement of Various Pulmonary Supplies	Qty	Unit	Brand	Statement of Compliance	
has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes single patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL full face mask, vented with anti asphyxiation valve has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes single patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a dual swivel elbow can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE full face mask, vented with anti asphyxiation valve has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes indicate the sum of the provided for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE full face mask, vented with anti asphyxiation valve has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes isingle patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow	1	BIPAP FULL FACE MASK (NON- VENTED)- MEDIUM	2	PCS			
nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		full face mask, vented with anti asphyxiation valve				t'-	
* material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow * can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL * full face mask, vented with anti asphyxiation valve * has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes * single patient use up to 14 days of use * material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow * can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE * full face mask, vented with anti asphyxiation valve * has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes * single patient use up to 14 days of use * material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow							
polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		single patient use up to 14 days of use				1	
* has a dual swivel elbow * can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL * full face mask, vented with anti asphyxiation valve * has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes * single patient use up to 14 days of use * material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow * can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE * full face mask, vented with anti asphyxiation valve * has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes * single patient use up to 14 days of use * material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow						1	
*can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL 3 PCS *full face mask, vented with anti asphyxiation valve *has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes *single patient use up to 14 days of use *material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam *has a simple color coded design to help choose the right mask size *has a dual swivel elbow *can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE 1 PCS *full face mask, vented with anti asphyxiation valve *has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes *single patient use up to 14 days of use *material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam *has a simple color coded design to help choose the right mask size *has a dual swivel elbow		has a simple color coded design to help choose the right mask size					
Note: Please provide 1 sample for evaluation 2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL 3 PCS - full face mask, vented with anti asphyxiation valve - has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes - single patient use up to 14 days of use - material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam - has a simple color coded design to help choose the right mask size - has a dual swivel elbow - can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE - full face mask, vented with anti asphyxiation valve - has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes - single patient use up to 14 days of use - material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam - has a simple color coded design to help choose the right mask size - has a dual swivel elbow		has a dual swivel elbow					
2 BIPAP FULL FACE MASK (NON- VENTED)- SMALL • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		can withstand pressures up to 25 cmH20					
• full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		Note: Please provide 1 sample for evaluation					
has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes single patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation BIPAP FULL FACE MASK (VENTED)- LARGE full face mask, vented with anti asphyxiation valve has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes single patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow	2	BIPAP FULL FACE MASK (NON- VENTED)- SMALL	3	PCS			
• has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		full face mask, vented with anti asphyxiation valve					
* material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow * can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE * full face mask, vented with anti asphyxiation valve * has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes * single patient use up to 14 days of use * material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam * has a simple color coded design to help choose the right mask size * has a dual swivel elbow							
polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		• single patient use up to 14 days of use					
• has a dual swivel elbow • can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow							
• can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		has a simple color coded design to help choose the right mask size					
Note: Please provide 1 sample for evaluation 3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		has a dual swivel elbow					
3 BIPAP FULL FACE MASK (VENTED)- LARGE • full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		can withstand pressures up to 25 cmH20					
• full face mask, vented with anti asphyxiation valve • has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		Note: Please provide 1 sample for evaluation					
has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes single patient use up to 14 days of use material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow	3	BIPAP FULL FACE MASK (VENTED)- LARGE	1	PCS			
nasogastric or orogastric tubes • single patient use up to 14 days of use • material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow		full face mask, vented with anti asphyxiation valve					
material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam has a simple color coded design to help choose the right mask size has a dual swivel elbow							
polyester, polyurethane foam • has a simple color coded design to help choose the right mask size • has a dual swivel elbow							
has a dual swivel elbow							
		has a simple color coded design to help choose the right mask size					
can withstand pressures up to 25 cmH20		has a dual swivel elbow					
		can withstand pressures up to 25 cmH20				134	
Note: Please provide 1 sample for evaluation		Note: Please provide 1 sample for evaluation					

Item No.	Specification Procurement of Various Pulmonary Supplies		Statement of Compliance [Bidders must state here either "Comply" or "Not Comply" against each of the individual parameters of each Specification stating the corresponding performance parameter of the equipment offered. Statements of "Comply" or "Not Comply" must be supported by evidence in a Bidders Bid and cross-referenced to that evidence. Evidence shall be in the form of manufacturer's un-amended sales literature, unconditional statements of specification and compliance issued by the manufacturer, samples, independent test data etc., as appropriate. A statement that is not supported by evidence or is subsequently found to be contradicted by the evidence presented will render the Bid under evaluation liable for rejection. A statement either in the Bidder's statement of compliance or the supporting evidence that is found to be false either during Bid evaluation, post-qualification or the execution of the Contract may be regarded as fraudulent and render the Bidder or supplier liable for prosecution subject to the applicable laws and issuances.]					
			Unit	Brand	Statement of Compliance			
4	BIPAP FULL FACE MASK (VENTED)- MEDIUM	38	PCS					
	full face mask, vented with anti asphyxiation valve		7.11					
	has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes							
	• single patient use up to 14 days of use							
	• material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam							
	has a simple color coded design to help choose the right mask size							
	has a dual swivel elbow							
	can withstand pressures up to 25 cmH20							
	Note: Please provide 1 sample for evaluation							
_								
5	BIPAP FULL FACE MASK (VENTED)- SMALL	38	PCS					
	full face mask, vented with anti asphyxiation valve							
	 has a tube fit design that maintan an effective seal around nasogastric or orogastric tubes 							
	 single patient use up to 14 days of use 							
	 material made up of silicone, acetal, polycarbonate nylon, polyester, polyurethane foam 							
	has a simple color coded design to help choose the right mask size			1				
	has a dual swivel elbow .							
	• can withstand pressures up to 25 cmH20 Note: Please provide 1 sample for evaluation			<u> </u>				
-V/29								
6	HEAT AND MOISTURE EXCHANGER (ADULT)	611	PCS	¥ ¥				
	 with mechanical Hydrophobic membrane 							
	with hygroscopic membrane							
	Individually packaged				30			
	 designed for tidal volume 150-1500ml 							
	Note: Please provide 1 sample for evaluation				BLY			
7	HEAT AND MOISTURE EXCHANGER (PEDIA)	300	PCS					
	 with mechanical Hydrophobic membrane 							
	with hygroscopic membrane							
	Individually packaged							
	 designed for tidal volume 75-500ml 	And Viller						
	Note: Please provide 1 sample for evaluation							

Item	ed as indicated in the Bid Bulletin No. 01-2023) Specification	Statement of Compliance						
No.			[Bidders must state here either "Comply" or "Not Comply" against each the individual parameters of each Specification stating the corresponding performance parameter of the equipment offered. Statements of "Comply" on "Not Comply" must be supported by evidence in a Bidders Bid and cross referenced to that evidence. Evidence shall be in the form of manufacturer unamended sales literature, unconditional statements of specification and compliance issued by the manufacturer, samples, independent test data etcas appropriate. A statement that is not supported by evidence or subsequently found to be contradicted by the evidence presented will render the Bid under evaluation liable for rejection. A statement either in the Bidder's statement of compliance or the supporting evidence that is found be false either during Bid evaluation, post-qualification or the execution the Contract may be regarded as fraudulent and render the Bidder or supplied liable for prosecution subject to the applicable laws and issuances.]					
	Procurement of Various Pulmonary Supplies	Qty	Unit	Brand	Statement of Compliance			
8	METERED DOSE INHALER PORT	127	PCS					
	• diameter- 22 mm (ID), 22 mm (OD)	- X						
	 with arrow that clearly marks the direction of aiflow 							
	disposable		^		the state of the s			
	 use in the delivery of aerosolized medication 	3	•					
	individually packed		y.	1 7 9 2				
	clear and transparent plastic							
	Note: Please provide 1 sample for evaluation							
9	PULSE OXIMETER SENSOR, FINGER PROBE (ADULT)	80	PCS					
	disposable							
	for single use							
	Intended for non-invasive, continuous SpO2 monitoring							
	individually packed							
	Ideal for application on the index finger, thumb, or great toe							
	Note: Please provide 1 sample for evaluation							
10	PULSE OXIMETER SENSOR, FINGER PROBE (NEONATE)	80	PCS		W P			
	disposable							
	for single use	y.						
	individually packed				1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1			
	 Intended for non-invasive, continuous SpO2 monitoring 				I Van			
	Ideal for application on the index finger, thumb, or great toe							
	Note: Please provide 1 sample for evaluation							
11	PULSE OXIMETER SENSOR, FINGER PROBE (PEDIA)	80	PCS					
	disposable							
	for single use							
	individually packed							
	Intended for non-invasive, continuous SpO2 monitoring							
	Ideal for application on the index finger, thumb, or great toe			****				
	Note: Please provide 1 sample for evaluation							
12	VIRAL FILTER	750	PCS					
	Individually packaged							
	designed for tidal volume 150-1200ml							
	• viral efficiency 99.99+%							
100	• 22 mm ID x 15 mm ID/22 mm OD							
	Bacterial efficiency 99.999+%			15				
	Note: Please provide 1 sample for evaluation							

Item No.	Specification	Statement of Compliance [Bidders must state here either "Comply" or "Not Comply" against each of the individual parameters of each Specification stating the corresponding performance parameter of the equipment offered. Statements of "Comply" of "Not Comply" must be supported by evidence in a Bidders Bid and cross referenced to that evidence. Evidence shall be in the form of manufacturer un-amended sales literature, unconditional statements of specification and compliance issued by the manufacturer, samples, independent test data etc as appropriate. A statement that is not supported by evidence or subsequently found to be contradicted by the evidence presented will render the Bid under evaluation liable for rejection. A statement either in the Bidder's statement of compliance or the supporting evidence that is found to be false either during Bid evaluation, post-qualification or the execution of the Contract may be regarded as fraudulent and render the Bidder or supplied liable for prosecution subject to the applicable laws and issuances.]					
	Procurement of Various Pulmonary Supplies	Qty	Unit	Brand	Statement of Compliance		
13	VOLUME RESPIRATOR TUBINGS (ADULT)	500	PCS				
	disposable						
	• for single use		1000				
	with nebulization kit						
	with water trap						
	• dual limb						
	Note: Please provide 1 sample for evaluation						
14	VOLUME RESPIRATOR TUBINGS (PEDIA)	250	PCS				
	disposable						
	• for single use						
	with nebulization kit						
	with water trap						
	• dual limb						
	Note: Please provide 1 sample for evaluation						
15	VOLUME RESPIRATOR TUBINGS (NEONATE)	300	PCS				
	disposable						
	• for single use				•		
	with nebulization kit						
	with water trap						
	dual limb			393			
	Note: Please provide 1 sample for evaluation						
	Name of Bidder: Submitted by: (Printed Name and Signature) Designation: Date:	-					